									Server					T.
Format	Extensions	Result	Big?	Folder in test images good	Folder in from_skyking	Notes	Date	Tester	(howe / gretzky)	Build Version•	OS if relevant	Summary of what was tested	Documents/Screenshots on squig**	Comment / Bug / Other?
AIM	aim		No		aim	110100			g , ,	- Lana Foreign		Canimary or innat true toolea	940.9	Commont Bug / Culor
						OME-XML file doesn't								
						contain any BinData,								
Al' ALOD	-10-1			-10-1	-10-1	so it is expected to be								
Alicona AL3D	al3d		No	al3d	al3d	black File: hiv.								
	am,					Expected to be mostly								
	amiramesh,					black, but there should								
	grey, hx,					be data between Z =								
Amira	labels		No		amira	140 and 149								
Analyze 7.5	img, hdr		No	analyze	analyze									
Andor	tif		No	andor										
						This is hard to view.								
						There are a couple of *really* bright pixels								
						that throw off the								
						display range; if you								
						set the min/max values								
						to roughly 600-800 (the	Э							
						real min/max for the								
Andor SIF	sif		No	andor-sif	andor-sif	first plane) then it should look better								
Ariol	jpeg		No	ariol	undor-sii	one file = a pyramid								
ARF	arf		No		samples/ARFTest	2.10 mo a pyramiu								
Audio Video Interleave			No		avi									
BD Pathway	exp, tif		Plate		bd-pathway									
Bio-Rad GEL	1sc		No	biorad-gel	biorad-gel									
				biorad-pic, biorad-										
Bio-Rad PIC	pic, xml, raw		No	raw	samples/ome, biorad									
Kodak Molecular Imaging	bip		No	bip	bip									
Windows Bitmap	bmp		No		bmp									
willdows bitiliap	ыпр		INO	ыпр	DITIP	picking any one of the								
						'acqp' files should								
						work, e.g.								
						bruker/wis/RA_060110								
						This is not a plate, but								
Bruker	acqp		No	bruker	bruker	rather an MRI acquisition.								
Burleigh	img		No		sxm/steve/Burleigh	acquisition.								
Duncign	cr2, crw, jpg,		110	buildigii	3XIII/3tcvc/Dullcigii									
Canon RAW	thm, wav		No	canon	canon									
	cr2, crw, jpg,													
	thm, wav, tif,													
DNG	tiff		No Plate		canon									
Cellomics C01	c01, dib apl, tnb, mtb,		Plate	cellomics	cellomics									
Olympus APL	tif		No	cellr	cellr									
						Takes ages. Uses								
CellSens VSI	vsi, ets		Yes		cellsens	sub-resolutions.								
CellWorx	pnl, htd, log		Plate	cellworx	cellworx									
	dic, dcm,													
	dicom, jp2, j2ki, j2kr, raw													
DICOM	ima	,	No	dicom	dicom									
D.00.	dv, r3d,		110	dicom	diodiii									
	r3d_d3d, dv.													
Deltavision	log, r3d.log		No		deltavision									
ECAT7	V		No	ecat7	ecat7									
Encapsulated PostScript	eps, epsi, ps		No	eps	eps									
FEI TIFF	tif, tiff		No		fei									
FEI/Philips	img		No		fei									
Flexible Image	9		T											
Transport System	fits, fts		No		fits									
Evotec Flex	flex, mea, res	3	Plate	flex	flex									
Olympus														
Fluoview/ABD TIFF	tif, tiff		No		fluoview, andor									
Fuji LAS 3000	img, inf oib, oif, pty,		No	fuji	fuji									
Olympus FV1000	lut		No	fv1000	fv1000									
Gatan Digital	1.4													
Micrograph	dm3		Usually - 4K x 4K	gatan	gatan									
Gatan DM2	dm2		Usually - 4K x 4K	gatan-dm2	gatan-dm2									
Amersham Biosciences GEL	gol		No	gol	anl									
Graphics Interchange	gel		INU	gel	gel									
Format	gif		No	gif	gif									
	1				-	I don't recommend								
						bothering with NDPI or								
						NDPIs, as these								
						formats are known to								
Hamamatsu NDPI	ndpi		Yes	hamamatsu	hamamatsu	be horribly limited. Uses sub-resolutions.								
			. 00			_ 500 005 100010110115.								
namamatsu NDF1			No		sxm/steve/Hamamatsu									

									0					
				Folder in					Server (howe /				Documents/Screenshots on	
Format	Extensions	Result	Big?	test images good	Folder in from_skyking	Notes	Date	Tester	gretzky)	Build Version•	OS if relevant	Summary of what was tested	squig**	Comment / Bug / Other?
1 Omiat	Exteriorie	rtoouit	Dig.	toot_imagoo_good	r older in from_okylang	File not valid, see	Duto	100001	growny	Build Voluloi1	GG II TOIGTUIN	Cultillary of what was tooled	oquig	Commont / Bug / Culor :
						https://trac.								
						openmicroscopy.org.								
Hamamatsu VMS	vms		Yes	hamamatsu-vms	hamamatsu-vms	uk/ome/ticket/9596								
						I don't recommend								
						bothering with NDPI or								
						NDPIs, as these formats are known to								
Hamamatsu NDPIS	ndpis		Yes	hamamatsu/philippe	hamamatsu/philippe	be horribly limited.								
Hitachi	txt			hitachi	hitachi	be normbly limited.								
Image Cytometry	LAL .		110	Tillacili	Intacin									
Standard	ics, ids		No	ics	ics									
Imacon	fff		No	imacon	imacon									
Image-Pro Sequence	seq		No	imagepro	imagepro									
Image-Pro Workspace				imagepro	imagepro									
IMAGIC	hed			imagic	imagic									
						Use the "file"								
						command to tell the								
						difference. Imaris 3								
						files will be reported as	5							
						"TIFF image data";								
						Imaris 5.5 files will be reported as								
						"Hierarchical Data								
						Format (version 5)								
						data"; the older raw								
						Imaris format files will								
Bitplane Imaris 3						be reported as simply								
(TIFF)	ims		No	imaris	imaris/kay	"data".	1							
						Use the "file" command to tell the								
						difference. Imaris 3								
						files will be reported as								
						"TIFF image data";	<b>'</b>							
						Imaris 5.5 files will be								
						reported as								
						"Hierarchical Data								
						Format (version 5)								
						data"; the older raw								
						Imaris format files will								
Bitplane Imaris 5.5					incode/abaia incode/incode	be reported as simply "data". Uses sub-								
(HDF)	ims		No	imaris	imaris/chris, imaris/jason, imaris/james, imaris/oliver									
(ПОГ)	IIIIS		INU	IIIIdiis	imans/james, imans/onver	Use the "file"								
						command to tell the								
						difference. Imaris 3								
						files will be reported as	3							
						"TIFF image data";								
						Imaris 5.5 files will be								
						reported as								
						"Hierarchical Data								
						Format (version 5) data"; the older raw								
						Imaris format files will								
						be reported as simply								
Bitplane Imaris	ims		No	imaris	imaris/generated	"data".								
					<b>y</b>	test image good:								
						image is black. That is								
						expected. Probably no								
						required to consider it								
						but contains ROIs. Useful for ROIs								
IMOD	mod		No	imod	imod	support/display								
Improvision TIFF	tif, tiff					auppoi vuispidy	-		-					
Improvision HFF	ul, ull		INU	improvision	improvision	KJH YFPPML 2000_A			-					
						101061 1: will								
						generate black								
	xdce, xml, tiff	,				thumbnails, raw data								
InCell 1000/2000	tif, xlog, im		Plate	incell	incell	missing								
			No. Melissa It is											
			technically a plate											
			but we lack all of											
			the metadata to be											
			able to put all of the											
			fields/wells together into a											
InCell 3000	frm		plate	incell3000	incell3000									
INR	inr				inr		-		-					
IPLab	ipl	-	No	iplab	iplab		-		-					
IVision	ipm				ivision		-		-					
JEOL	dat, img, par		No		sxm/steve/JEOL		-	-	-					
JPEG				jpeg, ariol	jpeg		-		-					
5. LO	jpg, jpeg, jpe jp2, j2k, jpf,		.10	JPO9, UNIO	JF~8		-		-					
JPEG-2000	jpz, jzk, jpi, jpx		Sometimes	jpeg2000	jpeg2000	Uses sub-resolutions.								
					sxm/steve/JPK Instr		1							
JPK Instruments	jpk		No	jpk	(Christian L)									

Format	Extensions	Result	Big?	Folder in test_images_good	Folder in from_skyking	Notes	Date	Tester	Server (howe / gretzky)	Build Version•	OS if relevant	Summary of what was tested	Documents/Screenshots on squig**	Comment / Bug / Other?
Khoros XV	xv		No	khoros	khoros				.,			·		
Li-Cor L2D	I2d, scn, tif		No	I2d	I2d									
Leica	lei, tif, tiff, raw	,	No	lei, lei_martin_sp2										
Leica Image File				leica-lif, leica-lif-										
Format	lif			martin	leica-lif									
Leica SCN	scn		yes	leica-scn	leica-scn	Uses sub-resolutions								
Leica TCS TIFF	tif, tiff, xml		No	leica-tcs	leica-tcs									
LEO	sxm, tif, tiff		No	leo	sxm/steve/LEO (Ewald B)									
		Cannot create				lifetime-extended.fli: float (display not								
LI-FLIM	fli	thumbnail	No	li-flim	li-flim	supported)								
Laboratory Imaging	lim		No	lim	lim									
Metamorph TIFF	tif, tiff		No	metamorph	metamorph									
MIAS	tif, tiff, txt		Plate	mias	mias									
				micromanager,										
Micro-Manager	tif, tiff, txt, xm	l .	No		micromanager									
MINC MRI	mnc		No	minc	minc									
Multiple Network			N1.											
Graphics	mng		No	mng	mng sxm/steve/Molecular									
Molecular Imaging	stp		No	molecular-imaging										
Medical Research	mrc, st, ali,		. 40	oiccular-imaying	maging									1
Council	map, rec		No	mrc	mrc									
Minolta MRW	mrw		Yes	mrw	mrw									
Hamamatsu														
Aquacosmos	naf		No	naf	naf									
NIfTI	nii, img, hdr		No	nifti	nifti									
						Nikon TIFFs have the Software tag set to "EZ-C1", i.e. running "tiffinfo foo.tiff" should								
Nikon Elements TIFF	tif, tiff		No	nikon	nikon	produce a line "Software: EZ-C1"								
						Nikon NEF TIFFs have Make/Model tags that contain "NIKON" and an EXIF IFD, which								
Nikon NEF	nef, tif, tiff		No	nikon	nikon	should be visible by running "tiffinfo foo.nef								
	tif. tiff					running unimo loo.nei								
Nikon TIFF			No	nikon nikon-jp2, nd2	nikon									
Nikon ND2 NOAA-HRD Gridded	nd2		No	nikon-jp2, na2	nikon-jp2, nd2									
Data Format			No	noaa-hrd	noaa-hrd									
NRRD	nrrd, nhdr		No	nrrd	nrrd									
	ome.tif, ome.			ome-tiff, ome-tif-										
OME-TIFF	tiff		No (test as Plate?)	multi	ome-tiff	tubhiswt.ome to be								
OME-XML	ome		No (test as Plate?)	ome-xml	ome-xml	removed embryo2.liff results in								
Openiab LIFF	liff		No	openlab	openiab	two images: - one has correct (3 channel 33z), -The second Image is expected to have SizeC = 1. It looks wrong, but is correct with respect to the metadata that is in the file.								
Openlab RAW	raw		No	openlab-raw	openlab-raw									
					<u> </u>	Image from Test image	•							
					sxm/steve/Oxford Instr	good: expected value								
Oxford Instruments	top		No	oxford-instruments	,	is size: 90x1								
PCX	рсх		No	рсх	рсх	No physical pixel size for DP1190.tif and a different physical pixel size PT1173A001.tif than for other formats.								
Perkin Elmer	bala to		Na	perkinelmer-	and design and the state of the	Expected: tiff converted and								
Densitometer	hdr, img		No	densitometer perkinelmer-	i e	metadata lost								
PerkinElmer Operetta			Plate	operetta	perkinelmer-operetta									
Portable Gray Map	pgm		No	pgm	OMERO_Test_Images/PG	N	-		-					-
PICT	pict, pct		No	pict		the animated PNG is png/Matrix.png, so if you imported everything in the 'png' folder then it's all set								
Animated PNG	png		No	png		roider then it's all set								-
POV-Ray	df3 tif, tiff, cfg,		No	povray	povray				-					
Prairie TIFF	xml		No	prairie	prairie									
Adobe Photoshop	psd		No	psd	psd									

				Folder in					Server (howe /				Documents/Screenshots on	
Format	Extensions	Result	Big?		Folder in from_skyking	Notes	Date	Tester	gretzky)	Build Version•	OS if relevant	Summary of what was tested	squig**	Comment / Bug / Other?
Pyramid TIFF	N/A		yes	pyramid	test_per_commit/pyramid	Uses sub-resolutions.								
Quesant AFM	afm		No		sxm/steve/Quesant Instr									
QuickTime	mov		No		quicktime									
RHK Technologies	sm2, sm3		No		sxm/steve/RHK Tech									
SBIG			No		sxm/steve/SBIG									
Olympus ScanR	dat, xml, tif		Plate	scanr	scanr	M/a accessment to come								
SPCImage Data	sdt		No	sdt	sdt	We currently use channels to store the 256 "bins". In the viewer we only display one channel at a time and use a lifetime slider (L:)								
			No		sxm/steve/Seiko Instr	Slider (L.)								
Seiko Compix Simple-PCI	xqd, xqf cxd		No											
SimplePCI TIFF	tif. tiff		No	simplepci simplepci	simplepci									
Olympus SIS TIFF	tif, tiff		No	sinpiepci	simplepci sis									
Olympus Slidebook	sld, spl		No	slidebook	slidebook									
SM Camera	οιυ, δμι		No	sm-camera	sxm/steve/SM-Camera		-							
SPIDER	spi		No	spider	spider		-							
Metamorph STK	stk, nd, tif, tiff		Plate		metamorph									
Aperio SVS	svs		Yes		svs	Uses sub-resolutions.								
Truevision Targa	tga		No		targa	OGCG GUD TCGGIULIONG.								
Text	txt, csv		No		text									
Adobe Photoshop	bit, cor		110	tont	toxt									
TIFF	tif, tiff		No	tiff	tiff									
Tagged Image File	tif, tiff, tf2, tf8					Uses sub-resolutions								
Format	btf		No	tiff	tiff	for JPEG2000 data.								
TillVision	vws, pst, inf		No	tillvision	tillvision									
TopoMetrix	tfr, ffr, zfr, zfp, 2fl		No	topometrix	sxm/steve/TopoMetrix									
Trestle	tif		Yes		trestle									
rresue	LII		165	uesue	trestie	uint32 (display not								
UBM	pr3		No	ubm	sxm/steve/UBM (Rosset)	supported)								
	ano, cfg, csv,				,	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,								
PerkinElmer	htm, rec, tim, zpo, tif		No	ultraview	ultraview									
reikilicillei	2p0, tii		INU		sxm/steve/Unisoku									
Unisoku STM	hdr, dat		No		(Takashi)									
Ombona Omi			110	umoona	(Tanaom)	float (display not								
Varian FDF	fdf		No	varian-mri	varian-mri	correctly supported)								
VG SAM	dti		No		sxm/steve/VG SAM									
Visitech XYS	xys, html		No	vg-sam visitech (xys only)	(Karen) visitech									
Volocity Library	xys, num		INU	visitecti (xys offiy)	VISILECTI									
Clipping	acff		No	volocity	volocity									
	mvd2, aisf,													
Volocity Library	aiix, dat, atsf		No	volocity	volocity									
					sxm/steve/WA Technology/WA Tech (Vladimir)/TOPSystem II, sxm/steve/WA Technology/WA Tech									
			No		(Adrian Jones)									
Zeiss CZI	czi		No	zeiss-czi	zeiss-czi									
Zeiss Laser-Scanning				zeiss-lsm, zeiss-										
Microscopy	Ism, mdb		No		zeiss-lsm									
Zeiss TIFF	tif		no	zeiss-tiff	zeiss-tiff									
Zeiss Vision Image (ZVI)	zvi		No	zeiss-zvi	zeiss-zvi									
Zip	zip		No		zip		-							
Zih	zιμ		INO	ZIP	ZIP									